

SCOPETECH^{DF}

HEALTHCARE Endoscope Cleaning System

KNIGHT^{HCC}
CFS TECHNOLOGIES
WE MAKE INSTRUMENT CLEANING EASY™

KEY BENEFITS

- Automated Scope Inflation
- Integrated Detergent Pump doses the right amount and tracks detergent conditions
- 2 Flushing / Rinsing Ports with Single or Dual Mode operation
- Cleaning Sequence Control
- Scope Inflation and Cleaning Process Reports (Windows compatible software and Excel reports)
- Sequence and Express modes
- Antimicrobial, IP-55 rated enclosure mounts on wall or counter top
- Electric pump eliminates technician muscle fatigue

4-IN-1 FLEXIBLE ENDOSCOPE CLEANING SYSTEM

The ScopeTech DF™ is a comprehensive flexible endoscope cleaning aid system that dramatically improves the speed of endoscope processing while providing step-by-step cleaning process management, cycle validation, and reports.

ScopeTech DF offers an automated, hands-free method for immersion leak test, 2 flushing ports with pulsating flush technology to remove bioburden, precision detergent dosing to achieve optimum cleaning conditions, and reports for management and surveyor audits.

An optional barcode scanner tracks scope serial number and technician I.D.s to ease documentation and reporting.



ScopeTech DF settings can be programmed on the user interface or by using the Windows compatible programming/reporting software. Settings are uploaded to the system via a USB flash drive and logged cycle data is retrieved the same way.



TECHNICAL DATA

Dimensions / Weight

9 3/4" H x 7 1/8" W x 8 3/4" D
24.76 cm x 18.09 cm x 22.22 cm

Cabinet Materials; Case Rating

#304 Stainless Steel, Powder Coated;
IP-55

Flow Rates

Flush Pump – .3 GPM (1.13 LPM)
Dosing Pump – 10 oz./min (29.6 ml/min.)

Flush Pump Suction

Maximum suction height = 24"/61 cm

Inflation Test Pressure

2.5-3 PSI (.17 - .27 bar)

Power Supply/Voltage

Wall Mount Type, IN: 100-240 Volts AC
0.8A 50-60 Hz, Out: 24 Volts DC 1.25A

Chemical Compatibility

Flush Pump, Flush Tubes, Hook Ups – Industry standard enzymatic and other detergents in diluted form. Do not use Isopropyl Alcohol with the flush pump.

Dosing Pump – Industry standard enzymatic and other detergents.

Unit Weight

13.3 lbs. (6.06 Kg)

Temperature Sensor

Max temperature = 160°C

Fuse – 2 Amp

Approvals

CSAus, CE Standard: CAN/CSA-C22.2
No. 61010-1 UL Std. No. 61010-1 IEC
61010-1

Product Testing

CSA International, Irvine CA

Part #	Mode
<i>ScopeTech DF models listed below include a flow meter, temperature sensor, low voltage power supply, international plug adapters, and accessory kit.</i>	
7260010	ScopeTech DF Flexible Endoscope Cleaning System with Barcode Scanner
7260011	ScopeTech DF Flexible Endoscope Cleaning System <i>without</i> barcode scanner
7117161-0	Inflation Test Kit, Olympus
0600818	Inflation Test Kit with Backflow Preventer, Pentax
7117161-F	Inflation Test Kit, Fujinon
7117161-P	Inflation Test Kit, Pentax
7117161-S	Inflation Test Kit, Karl Storz



7117161-0

0600818

7117161-F

7117161-P

7117161-S

ADDITIONAL DETAILS


Operating Modes

- Advanced "Sequence" mode requires the input of a personal I.D. and endoscope serial number to start the inflation test and cleaning operation. Technicians then step through the timed cleaning sequence following the programmed sequence.
- Basic "Express" mode is for pre-cleaning operations where speed or simplicity of the pre-cleaning operation are preferred or where cleaning of two endoscopes at once is required.

Cycle validation reports

On the User Interface

- ScopeTech DF reports can be viewed on the user display and provide managers fast and easy insight into scope pre-cleaning activities and performance. Records are sorted by endoscope type/model and serial number.



On your PC

- Printable report content saved in ScopeTech DF is retrieved using a USB flash drive then uploaded to the Scope Tech Windows compatible software tool installed on a laptop or desktop PC. Reports provide management validation of the inflation test, channel flushing and detergent dosing operations as well as confirmation of dates, times and technicians responsible for processing. Reports can be exported to Excel for printing or sharing with others. SettinAlerts when detergent level is low

Barcode Scan Functions

- Scan technician names and endoscope serial numbers quickly and easily.
- PC software automatically generates barcodes for technicians and scopes.
- Print barcode pages from PC software and post in a booklet or on the wall above the sink for convenient scanning.